



CVE-2007-6756

- MITRE
- NVD
- CVE.ORG
- JSON API
- Print: PDF

Summary

CVE	CVE-2007-6756
State	PUBLISHED
Assigner	mitre
Source Priority	CVE Program / NVD first with legacy fallback
Published	2014-08-12 16:55:04 UTC
Updated	2026-05-06 22:30:45 UTC
Description	ZOLL Defibrillator / Monitor M Series, E Series, and R Series have a default password for System Configuration mode, which allows an attacker to gain access to the device's configuration and potentially execute arbitrary code.

Risk And Classification

Primary CVSS: v2.0 4.9 from nvd@nist.gov

AV:L/AC:L/Au:N/C:N/I:N/A:C

Problem Types: CWE-255 | n/a

CVSS v2.0 Breakdown

Access Vector

Local

Access Complexity

Low

Authentication

None

Confidentiality

None

Integrity

None

Availability

Complete

AV:L/AC:L/Au:N/C:N/I:N/A:C

NVD Known Affected Configurations (CPE 2.3)

Type	Vendor	Product	Version	Update	Edition	Language
Hardware	Zoll	Monitor/defibrillator	e	All	All	All

Hardware	Zoll	Monitor/defibrillator	m	All	All	All
Hardware	Zoll	Monitor/defibrillator	r	All	All	All
Vendor Declared Affected Products						
Source	Vendor	Product	Version	Platforms		
CNA	Na	N/a	affected n/a	Not specified		
References						
Reference	Source		Link	Tags		
www.zoll.com/WorkArea/DownloadAsset.aspx	af854a3a-2127-422b-91ae-364da2661108		www.zoll.com			
www.zoll.com/WorkArea/DownloadAsset.aspx	af854a3a-2127-422b-91ae-364da2661108		www.zoll.com			
www.zoll.com/WorkArea/DownloadAsset.aspx	af854a3a-2127-422b-91ae-364da2661108		www.zoll.com			
IBM X-Force Exchange	af854a3a-2127-422b-91ae-364da2661108		exchange.xforce.ibmcloud.com			
CVE Program record	CVE.ORG		www.cve.org	canonical		
NVD vulnerability detail	NVD		nvd.nist.gov	canonical, and		
No vendor comments have been submitted for this CVE.						
There are currently no legacy QID mappings associated with this CVE.						