



# CVE-2018-7526

Published on: 05/24/2018 12:00:00 AM UTC

Last Modified on: 03/23/2021 11:24:38 PM UTC

## CVE-2018-7526

Source: Mitre

Source: NIST

CVE.ORG

Print: PDF



Certain versions of [Scroll Medical Air Systems](#) from [Beaconmedaes](#) contain the following vulnerability:

In TotalAlert Web Application in BeaconMedaes Scroll Medical Air Systems prior to v4107600010.23, by accessing a specific uniform resource locator (URL) on the webserver, a malicious user may be able to access information in the application without authenticating.

CVE-2018-7526 has been assigned by [ics-cert@hq.dhs.gov](mailto:ics-cert@hq.dhs.gov) to track the vulnerability - currently rated as **HIGH** severity.

CVSS3 Score: **7.5 - HIGH**

| Attack Vector    | Attack Complexity      | Privileges Required | User Interaction    |
|------------------|------------------------|---------------------|---------------------|
| <b>NETWORK</b>   | <b>LOW</b>             | <b>NONE</b>         | <b>NONE</b>         |
| Scope            | Confidentiality Impact | Integrity Impact    | Availability Impact |
| <b>UNCHANGED</b> | <b>HIGH</b>            | <b>NONE</b>         | <b>NONE</b>         |

CVSS2 Score: **5 - MEDIUM**

| Access Vector          | Access Complexity | Authentication      |
|------------------------|-------------------|---------------------|
| <b>NETWORK</b>         | <b>LOW</b>        | <b>NONE</b>         |
| Confidentiality Impact | Integrity Impact  | Availability Impact |
| <b>PARTIAL</b>         | <b>NONE</b>       | <b>NONE</b>         |

## CVE References

| Description                                               | Tags                                                                                                                                                                              | Link                                                                                                                       |
|-----------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| BeaconMedaes TotalAlert Scroll Medical Air Systems   CISA | <a href="#">Mitigation</a><br><a href="#">Third Party Advisory</a><br><a href="#">US Government Resource</a><br><a href="#">ics-cert.us-cert.gov</a><br><a href="#">text/html</a> | <a href="https://ics-cert.us-cert.gov/advisories/ICSMA-18-144-01">MISC ics-cert.us-cert.gov/advisories/ICSMA-18-144-01</a> |

There are currently no QIDs associated with this CVE

| Type                                                                                              | Vendor                       | Product                                             | Version | Update | Edition | Language |
|---------------------------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------|---------|--------|---------|----------|
| Hardware         | <a href="#">Beaconmedaes</a> | <a href="#">Scroll Medical Air Systems</a>          | -       | All    | All     | All      |
| Hardware         | <a href="#">Beaconmedaes</a> | <a href="#">Scroll Medical Air Systems</a>          | -       | All    | All     | All      |
| Operating System                                                                                  | <a href="#">Beaconmedaes</a> | <a href="#">Scroll Medical Air Systems Firmware</a> | All     | All    | All     | All      |
| Operating System                                                                                  | <a href="#">Beaconmedaes</a> | <a href="#">Scroll Medical Air Systems Firmware</a> | All     | All    | All     | All      |
| <pre>cpe:2.3:h:beaconmedaes:scroll_medical_air_systems:-:~::~~::~~::~~::~~::~~::~~::</pre>        |                              |                                                     |         |        |         |          |
| <pre>cpe:2.3:h:beaconmedaes:scroll_medical_air_systems:-:~::~~::~~::~~::~~::~~::~~::</pre>        |                              |                                                     |         |        |         |          |
| <pre>cpe:2.3:o:beaconmedaes:scroll_medical_air_systems_firmware:~::~~::~~::~~::~~::~~::~~::</pre> |                              |                                                     |         |        |         |          |
| <pre>cpe:2.3:o:beaconmedaes:scroll_medical_air_systems_firmware:~::~~::~~::~~::~~::~~::~~::</pre> |                              |                                                     |         |        |         |          |

[← Previous ID](#)

Next ID→

Use of this information constitutes acceptance for use in an AS IS condition. There are NO warranties, implied or otherwise, with regard to this information or its use. Any use of this information is at the user's risk. It is the responsibility of user to evaluate the accuracy, completeness or usefulness of any information, opinion, advice or other content. EACH USER WILL BE SOLELY RESPONSIBLE FOR ANY consequences of his or her direct or indirect use of this web site. ALL WARRANTIES OF ANY KIND ARE EXPRESSLY DISCLAIMED. This site will NOT BE LIABLE FOR ANY DIRECT, INDIRECT or any other kind of loss.

CVE.report and Source URL Uptime Status [status.cve.report](https://status.cve.report)