




CVE-2020-15485

Published on: 08/26/2020 12:00:00 AM UTC

Last Modified on: 07/15/2021 07:42:00 PM UTC

CVE-2020-15485

[Source: Mitre](#)[Source: NIST](#)[CVE.ORG](#)[Print: PDF](#) 

Certain versions of [M1000 Multipara Patient Monitor](#) from [Niscomed](#) contain the following vulnerability:

An issue was discovered on Nescomed Multipara Monitor M1000 devices. The onboard Flash memory stores data in cleartext, without integrity protection against tampering.

CVE-2020-15485 has been assigned by [M](#) cve@mitre.org to track the vulnerability - currently rated as **MEDIUM** severity.

CVSS3 Score: **5.5 - MEDIUM**

Attack Vector	Attack Complexity	Privileges Required	User Interaction
LOCAL	LOW	LOW	NONE
Scope	Confidentiality Impact	Integrity Impact	Availability Impact
UNCHANGED	HIGH	NONE	NONE

CVSS2 Score: **2.1 - LOW**

Access Vector	Access Complexity	Authentication
LOCAL	LOW	NONE
Confidentiality Impact	Integrity Impact	Availability Impact
PARTIAL	NONE	NONE

CVE References

Description	Tags	Link
Multipara Monitor - Multipara Monitor (Aqua 8) Manufacturer from New Delhi	Product Vendor Advisory www.niscomed.com text/html	AE MISC www.niscomed.com/multipara-monitor.html

By selecting these links, you may be leaving CVEreport webspace. We have provided these links to other websites because they may have information that would be of interest to you. No inferences should be drawn on account of other sites being referenced, or not, from this page. There may be other websites that are more appropriate for your purpose. CVEreport does not necessarily endorse the views expressed, or concur with the facts presented on these sites. Further, CVEreport does not endorse any commercial products that may be mentioned on these sites. Please address comments about any linked pages to comment@cve.report.

There are currently no QIDs associated with this CVE

Known Affected Configurations (CPE V2.3)

Type	Vendor	Product	Version	Update	Edition	Language
Hardware 	Niscomed	M1000 Multipara Patient Monitor	-	All	All	All
Operating System	Niscomed	M1000 Multipara Patient Monitor Firmware	-	All	All	All
Hardware 	Niscomed	Multipara Monitor M1000	-	All	All	All
Hardware 	Niscomed	Multipara Monitor M1000	-	All	All	All
Operating System	Niscomed	Multipara Monitor M1000 Firmware	-	All	All	All
Operating System	Niscomed	Multipara Monitor M1000 Firmware	-	All	All	All
cpe:2.3:h:niscomed:m1000_multipara_patient_monitor:-:*:*:*:*:*:						
cpe:2.3:o:niscomed:m1000_multipara_patient_monitor_firmware:-:*:*:*:*:*:						
cpe:2.3:h:niscomed:multipara_monitor_m1000:-:*:*:*:*:*:						
cpe:2.3:h:niscomed:multipara_monitor_m1000:-:*:*:*:*:*:						
cpe:2.3:o:niscomed:multipara_monitor_m1000_firmware:-:*:*:*:*:*:						
cpe:2.3:o:niscomed:multipara_monitor_m1000_firmware:-:*:*:*:*:*:						

No vendor comments have been submitted for this CVE

Social Mentions

Source	Title	Posted (UTC)
← Previous IDNext ID →		

site will NOT BE LIABLE FOR ANY DIRECT, INDIRECT or any other kind of loss.

CVE, CWE, and OVAL are registered trademarks of [The MITRE Corporation](#) and the authoritative source of CVE content is [MITRE's CVE web site](#). This site includes MITRE data granted under the following [license](#).

CVE.report and Source URL Uptime Status [status.cve.report](#)