




CVE-2020-27260

Published on: 01/08/2021 12:00:00 AM UTC

Last Modified on: 03/23/2021 11:23:28 PM UTC

CVE-2020-27260

[Source: Mitre](#)[Source: NIST](#)[CVE.ORG](#)[Print: PDF](#) 

Certain versions of [Vital Signs Monitor Vc150](#) from [Innokasmedical](#) contain the following vulnerability:

Innokas Yhtymä Oy Vital Signs Monitor VC150 prior to Version 1.7.15 HL7 v2.x injection vulnerabilities exist in the affected products that allow physically proximate attackers with a connected barcode reader to inject HL7 v2.x segments into specific HL7 v2.x messages via multiple expected parameters.

CVE-2020-27260 has been assigned by ics-cert@hq.dhs.gov to track the vulnerability - currently rated as **MEDIUM** severity.

CVSS3 Score: **5.3 - MEDIUM**

Attack Vector	Attack Complexity	Privileges Required	User Interaction
PHYSICAL	LOW	NONE	NONE
Scope	Confidentiality Impact	Integrity Impact	Availability Impact
CHANGED	NONE	HIGH	NONE

CVSS2 Score: **2.1 - LOW**

Access Vector	Access Complexity	Authentication
LOCAL	LOW	NONE
Confidentiality Impact	Integrity Impact	Availability Impact
NONE	PARTIAL	NONE

CVE References

Description	Tags	Link
Innokas Yhtymä Oy Vital Signs Monitor CISA	Third Party Advisory US Government Resource us-cert.cisa.gov	MISC us-cert.cisa.gov/ics/advisories/icsma-21-007-01

By selecting these links, you may be leaving CVEreport webspace. We have provided these links to other websites because they may have information that would be of interest to you. No inferences should be drawn on account of other sites being referenced, or not, from this page. There may be other websites that are more appropriate for your purpose. CVEreport does not necessarily endorse the views expressed, or concur with the facts presented on these sites. Further, CVEreport does not endorse any commercial products that may be mentioned on these sites. Please address comments about any linked pages to comment@cve.report.

There are currently no QIDs associated with this CVE

Known Affected Configurations (CPE V2.3)

Type	Vendor	Product	Version	Update	Edition	Language
Hardware 	Innokasmedical	Vital Signs Monitor Vc150	-	All	All	All
Hardware 	Innokasmedical	Vital Signs Monitor Vc150	-	All	All	All
Operating System	Innokasmedical	Vital Signs Monitor Vc150 Firmware	All	All	All	All
Operating System	Innokasmedical	Vital Signs Monitor Vc150 Firmware	All	All	All	All
cpe:2.3:h:innokasmedical:vital_signs_monitor_vc150:-:*:*:*:*:*:						
cpe:2.3:h:innokasmedical:vital_signs_monitor_vc150:-:*:*:*:*:*:						
cpe:2.3:o:innokasmedical:vital_signs_monitor_vc150_firmware:*:*:*:*:*:						
cpe:2.3:o:innokasmedical:vital_signs_monitor_vc150_firmware:*:*:*:*:*:						

No vendor comments have been submitted for this CVE

[← Previous ID](#)

Next ID→